

Nur für Forschungszwecke

PARK2/Parkin Monoklonaler Antikörper

Katalog-Nr.: CL594-66674



Allgemeine Informationen

Katalog-Nr.:
CL594-66674

Größe:
100ul

Wirt:
Maus

Isotyp:
IgG2b

Immunogen Katalognummer:
AG5179

GenBank-Zugangsnummer:
BC022014

GeneID (NCBI):
5071

Vollständiger Name:
Parkinson disease (autosomal
recessive, juvenile) 2, parkin

Berechnete Masse:
52 kDa

Beobachtete Masse:
42-52 kDa

Reinigungsmethode:
Protein-A-Reinigung

CloneNo.:
2H5A7

Anregungs-/Emissionsmaxima-
Wellenlängen:
588 nm / 604 nm

Anwendungen

Geprüfte Anwendungen:

Getestete Reaktivität:
Human, Maus

Hintergrundinformationen

Parkin, a RING-type E3 ubiquitin-protein ligase, is involved in the ubiquitination pathway and contributes to protection from neurotoxicity induced by unfolded protein stresses. Its ubiquitin-protein ligase activity promotes the degradation of a variety of proteins including itself. Mutations in Parkin are implicated in the pathogenesis of autosomal recessive familial Parkinson's disease. It has 8 isoforms produced by alternative splicing.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

