

Nur für Forschungszwecke

Phospho-RIPK1 (Ser161) Monoklonaler Antikörper

Katalog-Nr.:CL594-66854



Allgemeine Informationen

Katalog-Nr.: CL594-66854	GenBank-Zugangsnummer: NM_003804	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von8737	GeneID (NCBI): von8737	CloneNo.: 1B2G1
Nanodrop;	Vollständiger Name: receptor (TNFRSF)-interacting serine- threonine kinase 1	Anregungs-/Emissionsmaxima- Wellenlängen: 588 nm / 604 nm
Wirt: Maus	Berechnete Masse: 76 kDa	
Isotyp: IgG1	Beobachtete Masse: 74 kDa, 45 kDa	

Anwendungen

Geprüfte Anwendungen:
FC (Intra)

Getestete Reaktivität:
Human

Hintergrundinformationen

RIPK1(Receptor-interacting serine/threonine-protein kinase 1) is primarily involved in mediating TNF-R1-induced cell activation, apoptosis and necroptosis and belongs to a novel class of kinases that function in cell survival and cell death mechanisms(PMID:22685397).It has 2 isoforms produced by alternative splicing.It also can exist as the protein with lower molecular weight about 25 kDa and 45 kDa in HAEC and HUVEC (PMID: 22685397).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.

Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ HEK-293T cells untreated (dashed lines) or treated with Calyculin A (red) were intracellularly stained with 0.13 ug CoraLite® 594 Anti-Human Phospho-RIPK1 (Ser161) (CL594-66854, Clone:1B2G1) (red), or 0.13 ug Control Antibody (blue). Cells were fixed with 4% PFA and permeabilized with 90% MeOH.