

Nur für Forschungszwecke

PDS5A Monoklonaler Antikörper

Katalog-Nr.:CL594-67520



Allgemeine Informationen

Katalog-Nr.:
CL594-67520

Größe:
100ul , Konzentration: 1000 µg/ml von23244
Nanodrop;

Wirt:
Maus

Isotyp:
IgG1

Immunogen Katalognummer:
AG11525

GenBank-Zugangsnummer:
BC009650

GeneID (NCBI):
1297

Vollständiger Name:
PDS5, regulator of cohesion
maintenance, homolog A (S.
cerevisiae)

Berechnete Masse:
1297 aa, 147 kDa

Reinigungsmethode:
Protein-G-Reinigung

CloneNo.:
1C8D4

Empfohlene Verdünnungen:
IF 1:50-1:500

Anregungs-/Emissionsmaxima-
Wellenlängen:
588 nm / 604 nm

Anwendungen

Geprüfte Anwendungen:
IF

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

IF : SMMC-7721-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed SMMC-7721 cells using CoraLite®594 PDS5A antibody (CL594-67520, Clone: 1C8D4) at dilution of 1:200.