

Nur für Forschungszwecke

CD11c Monoklonaler Antikörper

Katalog-Nr.:CL647-65086



Allgemeine Informationen

Katalog-Nr.: CL647-65086	GenBank-Zugangsnummer: BC038237	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 µl/test	GeneID (NCBI): 3687	CloneNo.: 3,9
Wirt: Maus	Vollständiger Name: integrin, alpha X (complement component 3 receptor 4 subunit)	Anregungs-/Emissionsmaxima-Wellenlängen: 654 nm / 674 nm
Isotyp: IgG1, kappa	Berechnete Masse: 1169 aa, 129 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11c, also known as integrin αX , is a type I transmembrane glycoprotein present on a variety of cells, including monocytes/macrophages, granulocytes, a subset of B cells, NK cells and dendritic cells (PMID: 2897326; 1680915; 1694698; 17389580). As a result of its high level of expression on most dendritic cells, CD11c is typically considered to be a marker of conventional dendritic cells (PMID: 27119555). CD11c forms an α/β heterodimer with CD18 (integrin $\beta 2$). CD11c/CD18 acts a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis (PMID: 1671533).

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



100 μ l human peripheral blood were surface stained with 5.00 μ l CoraLite® Plus 647-conjugated Anti-Human CD11c (CL647-65086, Clone:3.9) and then treated with red blood cell lysis buffer. Cells were not fixed.



100 μ l human peripheral blood were surface stained with CoraLite® Plus 647-conjugated Mouse IgG1 isotype control and then treated with red blood cell lysis buffer. Cells were not fixed.