

Nur für Forschungszwecke

# CD115 (c-fms) Monoklonaler Antikörper

Katalog-Nr.:CL647-65095



## Allgemeine Informationen

<b>Katalog-Nr.:</b> CL647-65095	<b>GenBank-Zugangsnummer:</b> BC036343	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100tests , 5 µl/test	<b>GeneID (NCBI):</b> 12978	<b>CloneNo.:</b> AFS98
<b>Wirt:</b> Ratte	<b>Vollständiger Name:</b> colony stimulating factor 1 receptor	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 654 nm / 674 nm
<b>Isotyp:</b> Rat IgG2a, kappa		

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Maus

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Bei 2-8°C lagern. Vor Licht schützen.

**Lagerungspuffer:**  
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> Balb/c mouse peritoneal macrophages were surface stained with 5 ul CoraLite® Plus 647-conjugated Anti-Mouse CD115 (c-fms) (CL647-65095, Clone: AF598) (blue) or CoraLite® Plus 647-conjugated Rat IgG2a isotype control (black). Cells were not fixed.