

Nur für Forschungszwecke

IL-12/IL-23 p40 Monoklonaler Antikörper

Katalog-Nr.:CL647-65097



Allgemeine Informationen

Katalog-Nr.:
CL647-65097

Größe:
100tests , 5 µl/test

Wirt:
Ratte

Isotyp:
Rat IgG2a, kappa

GenBank-Zugangsnummer:
BC103608

GeneID (NCBI):
16160

Vollständiger Name:
interleukin 12b

Reinigungsmethode:
Affinitätsreinigung

CloneNo.:
C17.8

Anregungs-/Emissionsmaxima-Wellenlängen:
654 nm / 674 nm

Anwendungen

Geprüfte Anwendungen:
FC (Intra)

Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



IFN- γ primed Balb/c mouse peritoneal macrophages stimulated with LPS in the presence of a protein transport inhibitor were treated with fixation and permeabilization reagents and intracellularly stained with CoraLite® Plus 647-Rat IgG2a Isotype Control (left) or 5 μ l CoraLite® Plus 647 Anti-Mouse IL-12/IL-23 p40 (CL647-65097, Clone:C17.8) (right).