

Nur für Forschungszwecke

CD24 Monoklonaler Antikörper

Katalog-Nr.:CL647-65121



Allgemeine Informationen

Katalog-Nr.:

CL647-65121

Größe:

100tests , 5 µl/test

Wirt:

Ratte

Isotyp:

IgG2b, kappa

GenBank-Zugangsnummer:

BC075622

GeneID (NCBI):

12484

Vollständiger Name:

CD24a antigen

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

M1/69

Anregungs-/Emissionsmaxima-
Wellenlängen:

654 nm / 674 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ C57BL/6 mouse splenocytes were surface stained with FITC-Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and 5.00 ul CoraLite® Plus 647-conjugated Anti-Mouse CD24 (CL647-65121, Clone: M1/69). Cells were not fixed.



1X10⁶ C57BL/6 mouse splenocytes were surface stained with FITC-Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and CoraLite® Plus 647-conjugated rat IgG2b isotype control antibody. Cells were not fixed.