

Nur für Forschungszwecke

CD4 Monoklonaler Antikörper

Katalog-Nr.:CL647-65199



Allgemeine Informationen

Katalog-Nr.:

CL647-65199

Größe:

100ug , 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2b, kappa

GenBank-Zugangsnummer:

BC039137

GeneID (NCBI):

12504

Vollständiger Name:

CD4 antigen

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

RM4-4

Anregungs-/Emissionsmaxima-
Wellenlängen:

654 nm / 674 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with CoraLite® 488 Plus Anti-Mouse CD3 (CL488-65077, Clone: 17A2) and 1 ug CoraLite® Plus 647 Anti-Mouse CD4 (CL647-65199, Clone: RM4-4) or CoraLite® Plus 647-Conjugated rat IgG2b isotype control. Cells were not fixed.