

Nur für Forschungszwecke

IgG1 Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL647-65212

Allgemeine Informationen

Katalog-Nr.:
CL647-65212
Größe:
100ug, 0.5 mg/mL
Wirt:
Ratte
Isotyp:
IgG1, κ

GenBank-Zugangsnummer:
GeneID (NCBI):
Vollständiger Name:

Reinigungsmethode:
Affinitätsreinigung
CloneNo.:
HRPN
Anregungs-/Emissionsmaxima-Wellenlängen:
654 nm / 674 nm

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:

Hintergrundinformationen

The HRPN immunoglobulin is useful as an isotype-matched control. The HRPN immunoglobulin has an unknown binding specificity and is used as an isotype control for rat IgG1 antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern. Vor Licht schützen.
Lagerungspuffer:
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 5 ul CoraLite® 488-conjugated Anti-Mouse CD4 (CL488-65104, Clone: GK1.5) and 5 ul CoraLite® Plus 647-conjugated Anti-Mouse CD25 (CL647-65137, Clone: PC61.5). Cells were not fixed.



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 5 ul CoraLite® 488-conjugated Anti-Mouse CD4 (CL488-65104, Clone: GK1.5) and CoraLite® Plus 647-conjugated Rat IgG1 Isotype Control (CL647-65212, Clone: HRPN). Cells were not fixed.