

Nur für Forschungszwecke

PLVAP Monoklonaler Antikörper

Katalog-Nr.:CL647-65214



Allgemeine Informationen

Katalog-Nr.:

CL647-65214

Größe:

100ug , 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC013517

GeneID (NCBI):

84094

Vollständiger Name:

plasmalemma vesicle associated
protein

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

MECA-32

Anregungs-/Emissionsmaxima-

Wellenlängen:

654 nm / 674 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ bEnd.3 cells were surface stained with 0.5 ug CoraLite® Plus 647 Anti-Mouse PLVAP (CL647-65214, Clone: MECA-32) (red) or 0.5 ug Control Antibody. Cells were not fixed.