

Nur für Forschungszwecke

# CD71 Monoklonaler Antikörper

Katalog-Nr.:CL647-65235



## Allgemeine Informationen

<b>Katalog-Nr.:</b> CL647-65235	<b>GenBank-Zugangsnummer:</b> BC001188	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100tests , 5 µl/test	<b>GeneID (NCBI):</b> 7037	<b>CloneNo.:</b> OKT9
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> transferrin receptor (p90, CD71)	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 654 nm / 674 nm
<b>Isotyp:</b> IgG1, kappa	<b>Berechnete Masse:</b> 85 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Lagerungspuffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



$1 \times 10^6$  PHA-activated human PBMCs were surface stained with 5  $\mu$ l CoraLite® Plus 647 Anti-Human CD71 (CL647-65235, Clone: OKT9), or APC Mouse IgG1 Isotype Control. Cells were not fixed.