

Nur für Forschungszwecke

# CD43 Monoklonaler Antikörper

Katalog-Nr.:CL647-65237



## Allgemeine Informationen

**Katalog-Nr.:**  
CL647-65237  
**Größe:**  
100ug , 500 µg/ml  
**Wirt:**  
Maus  
**Isotyp:**  
IgG1, kappa

**GenBank-Zugangsnummer:**  
**GeneID (NCBI):**  
24796  
**Vollständiger Name:**  
sialophorin  
**Berechnete Masse:**  
38 kDa

**Reinigungsmethode:**  
Affinitätsreinigung  
**CloneNo.:**  
W3/13  
**Anregungs-/Emissionsmaxima-Wellenlängen:**  
654 nm / 674 nm

## Anwendungen

**Geprüfte Anwendungen:**  
FC  
**Getestete Reaktivität:**  
Ratte

## Hintergrundinformationen


### Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
**Lagerungspuffer:**  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



$1 \times 10^6$  Wistar rat splenocytes were surface stained with 0.5  $\mu$ g CoraLite® Plus 647 Anti-Rat CD43 (CL647-65237, Clone: W3/13) or 0.5  $\mu$ g Control Antibody. Cells were not fixed.