

Nur für Forschungszwecke

# CD49b Monoklonaler Antikörper

Katalog-Nr.:CL647-65245



## Allgemeine Informationen

Katalog-Nr.:

CL647-65245

Größe:

100ug , 500 µg/ml

Wirt:

Armenian Hamster

Isotyp:

IgG

GenBank-Zugangsnummer:

BC065139

GeneID (NCBI):

16398

Vollständiger Name:

integrin alpha 2

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

HMA2

Anregungs-/Emissionsmaxima-

Wellenlängen:

654 nm / 674 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> mouse splenocytes were surface stained with FITC Anti-Mouse NK1.1 (CD161) (FITC-65138, Clone: PK136) and 0.5 ug CoraLite® Plus 647 Anti-Mouse CD49b (CL647-65245, Clone: HMA 2) or 0.5 ug CoraLite® Plus 647-conjugated Armenian Hamster IgG. Cells were not fixed.