

Nur für Forschungszwecke

# CD58 Monoklonaler Antikörper

Katalog-Nr.:CL750-65193



## Allgemeine Informationen

Katalog-Nr.:

CL750-65193

Größe:

100tests , 5 µl/test

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC005930

GeneID (NCBI):

965

Vollständiger Name:

CD58 molecule

Berechnete Masse:

28 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

TS2/9

Anregungs-/Emissionsmaxima-Wellenlängen:

755 nm / 780 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD58 (CL750-65193, Clone:TS2/9) or unstained. Cells were not fixed. Monocytes were gated.