

Nur für Forschungszwecke

CD7 Monoklonaler Antikörper

Katalog-Nr.:CL750-65203



Allgemeine Informationen

Katalog-Nr.: CL750-65203	GenBank-Zugangsnummer: BC009293	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 µl/test	GeneID (NCBI): 924	CloneNo.: 4H9
Wirt: Maus	Vollständiger Name: CD7 molecule	Anregungs-/Emissionsmaxima-Wellenlängen: 755 nm / 780 nm
Isotyp: Mouse IgG2a, kappa	Berechnete Masse: 240 aa, 25 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface co-stained with APC Anti-Human CD4 and 5 ul CoraLite® Plus 750 Anti-Human CD7 (CL750-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD7 (CL750-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.