

Nur für Forschungszwecke

IgG2b Isotype Control Monoklonaler Antikörper



Katalog-Nr.: CR-65128

Allgemeine Informationen

Katalog-Nr.:

CR-65128

Größe:

100ug, 0.2 mg/mL

Wirt:

Maus

Isotyp:

IgG2b

GenBank-Zugangsnummer:

GeneID (NCBI):

16016

Vollständiger Name:

immunoglobulin heavy chain 3
(serum IgG2b)

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

MPC-11

Anregungs-/Emissionsmaxima-
Wellenlängen:

592 nm / 611 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Hintergrundinformationen

The MPC-11 immunoglobulin is useful as an isotype-matched control. The MPC-11 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2b antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 0.2 ug Cardinal Red™ Anti-Mouse IgG2b Isotype Control (CR-65128, Clone:MPC-11) (red), 5 ul CR Anti-Human CD4 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.