

Nur für Forschungszwecke

CD4 Monoklonaler Antikörper

Katalog-Nr.: CR-65134



Allgemeine Informationen

Katalog-Nr.: CR-65134	GenBank-Zugangsnummer: BC025782	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests, 5 µl/test	GeneID (NCBI): 920	CloneNo.: OKT4
Wirt: Maus	Vollständiger Name: CD4 molecule	Anregungs-/Emissionsmaxima-Wellenlängen: 592 nm / 611 nm
Isotyp: IgG2b, kappa	Berechnete Masse: 55 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface co-stained with Atlantic Blue Anti-Human CD3 and 5 ul Cardinal Red™ Anti-Human CD4 (CR-65134, Clone:OKT4) or Mouse IgG2b Isotype Control. Cells were not fixed. Lymphocytes were gated.