

Nur für Forschungszwecke

CD7 Monoklonaler Antikörper

Katalog-Nr.: CR-65203



Allgemeine Informationen

Katalog-Nr.:

CR-65203

Größe:

100 tests, 5 µl/test

Wirt:

Maus

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC009293

GeneID (NCBI):

924

Vollständiger Name:

CD7 molecule

Berechnete Masse:

240 aa, 25 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

4H9

Anregungs-/Emissionsmaxima-

Wellenlängen:

592 nm / 611 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface co-stained with APC Anti-Human CD4 and 5 μ l Cardinal Red™ Anti-Human CD7 (CR-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 μ l Cardinal Red™ Anti-Human CD7 (CR-65203, Clone:4H9) (red) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.