

Nur für Forschungszwecke

CD95 Monoklonaler Antikörper

Katalog-Nr.: FITC-65071



Allgemeine Informationen

Katalog-Nr.:

FITC-65071

Größe:

100 tests, 10 µl/test

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC012479

GeneID (NCBI):

355

Vollständiger Name:

Fas (TNF receptor superfamily, member 6)

Berechnete Masse:

35-38 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

DX2

Anregungs-/Emissionsmaxima-

Wellenlängen:

494 nm / 520 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Fas (CD95/APO-1) is a transmembrane glycoprotein belonging to the tumor necrosis factor (TNF) receptor superfamily. It can mediate apoptosis by ligation with an agonistic anti-Fas antibody or Fas ligand. Stimulation of Fas results in the aggregation of its intracellular death domains, leading to the formation of the death-inducing signaling complex (DISC). FAS-mediated apoptosis may have a role in the induction of peripheral tolerance, in the antigen-stimulated suicide of mature T-cells, or both.

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood lymphocytes were surface stained with 10 ul FITC Anti-Human CD95 (FITC-65071, Clone:DX2) (green) or FITC-mouse IgG1 isotype control (black). Cells were not fixed.