

Nur für Forschungszwecke

CD115 (c-fms) Monoklonaler Antikörper

Katalog-Nr.: FITC-65095



Allgemeine Informationen

Katalog-Nr.:
FITC-65095

Größe:
100ug, 0.5 mg/ml

Wirt:
Ratte

Isotyp:
IgG2a, kappa

GenBank-Zugangsnummer:
BC036343

GeneID (NCBI):
12978

Vollständiger Name:
colony stimulating factor 1 receptor

Reinigungsmethode:
Affinitätsreinigung

CloneNo.:
AFS98

Anregungs-/Emissionsmaxima-
Wellenlängen:
494 nm / 520 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug FITC Anti-Mouse CD115 (c-fms) (FITC-65095, Clone: AFS98) (green) or 0.5 ug FITC-Rat IgG2a Isotype Control (black). Cells were not fixed.