

Nur für Forschungszwecke

CD1a Monoklonaler Antikörper

Katalog-Nr.: FITC-65172



Allgemeine Informationen

Katalog-Nr.:
FITC-65172

Größe:
100tests, 5 µl/test

Wirt:
Maus

Isotyp:
IgG1, kappa

GenBank-Zugangsnummer:
BC031645

GeneID (NCBI):
909

Vollständiger Name:
CD1a molecule

Berechnete Masse:
327 aa, 37 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit Fluoresceinisothiocyanat (FITC) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:
HI149

Anregungs-/Emissionsmaxima-Wellenlängen:
495 nm / 524 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1×10^6 MOLT-4 cells were surface stained with 5 μ l FITC Plus Anti-Human CD1a (FITC-65172, Clone:HI149) or Mouse IgG1 Isotype Control. Cells were not fixed.