

Nur für Forschungszwecke

CD58 Monoklonaler Antikörper

Katalog-Nr.: FITC-65193



Allgemeine Informationen

Katalog-Nr.:
FITC-65193

Größe:
100tests, 5 µl/test

Wirt:
Maus

Isotyp:
IgG1, kappa

GenBank-Zugangsnummer:
BC005930

GeneID (NCBI):
965

Vollständiger Name:
CD58 molecule

Berechnete Masse:
28 kDa

Reinigungsmethode:
Protein-G-Reinigung

CloneNo.:
TS2/9

Anregungs-/Emissionsmaxima-
Wellenlängen:
495 nm / 524 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 5 ul FITC Anti-Human CD58 (FITC-65193, Clone:TS2/9) or unstained. Cells were not fixed. Monocytes were gated.