

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: FITC-65208

Allgemeine Informationen

Katalog-Nr.:

FITC-65208

Größe:

100ug, 0.5 mg/mL

Wirt:

Maus

Isotyp:

IgG2a

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

C1.18.4

Anregungs-/Emissionsmaxima-Wellenlängen:

494 nm / 520 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Hintergrundinformationen

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ SJL mouse splenocytes were surface stained with 0.20 ug FITC Mouse IgG2a Isotype Control (C1.18.4) (FITC-65208, Clone: C1.18.4) (red) or with 0.20 ug FITC Anti-Mouse CD45.1 (FITC-65092, Clone: A20) (green), or not stained (blank) (blue). Cells were not fixed.