

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: FITC-65209

Allgemeine Informationen

Katalog-Nr.: FITC-65209	GenBank-Zugangsnummer: GeneID (NCBI):	Reinigungsmethode: Affinitätsreinigung
Größe: 100ug, 0.5 mg/mL	Vollständiger Name:	CloneNo.: 2A3
Wirt: Ratte		Anregungs-/Emissionsmaxima-Wellenlängen: 494 nm / 520 nm
Isotyp: IgG2a, kappa		

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:

Hintergrundinformationen

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.5 ug FITC Anti-Rat IgG2a Isotype Control (FITC-65209, Clone:2A3) (red), 0.5 ug FITC-65141 (green), or unstained (blue). Cells were not fixed.