

Nur für Forschungszwecke

# PLVAP Monoklonaler Antikörper

Katalog-Nr.: FITC-65214



## Allgemeine Informationen

Katalog-Nr.:

FITC-65214

Größe:

100ug , 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC013517

GeneID (NCBI):

84094

Vollständiger Name:

plasmalemma vesicle associated protein

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

MECA-32

Anregungs-/Emissionsmaxima-

Wellenlängen:

495 nm / 524 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> bEnd.3 cells were surface stained with 1 ug FITC Plus Anti-Mouse PLVAP (FITC-65214, Clone: MECA-32) (red) or 1 ug Control Antibody. Cells were not fixed.