

Nur für Forschungszwecke

# CD43 Monoklonaler Antikörper

Katalog-Nr.: FITC-65237



## Allgemeine Informationen

Katalog-Nr.:  
FITC-65237

Größe:  
100ug, 500 µg/ml

Wirt:  
Maus

Isotyp:  
IgG1, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):  
24796

Vollständiger Name:  
sialophorin

Berechnete Masse:  
38 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

W3/13

Anregungs-/Emissionsmaxima-  
Wellenlängen:

495 nm / 524 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Ratte

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



$1 \times 10^6$  Wistar rat splenocytes were surface stained with 0.5  $\mu$ g FITC Plus Anti-Rat CD43 (FITC-65237, Clone: W3/13) or 0.5  $\mu$ g Control Antibody. Cells were not fixed.