

Nur für Forschungszwecke

# CD94 Monoklonaler Antikörper

Katalog-Nr.:FITC-65240



## Allgemeine Informationen

Katalog-Nr.:

FITC-65240

Größe:

100ug , 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC120854

GeneID (NCBI):

16643

Vollständiger Name:

killer cell lectin-like receptor,  
subfamily D, member 1

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

18D3

Anregungs-/Emissionsmaxima-

Wellenlängen:

495 nm / 524 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> mouse splenocytes were surface co-stained with CoraLite®647 Anti-Mouse CD49b (HMA2) (CL647-65245, Clone: HMA2) and 0.5 µg FITC Plus Anti-Mouse CD94 (FITC-65240, Clone: 18D3) or 0.5 µg Isotype Control. Cells were not fixed.