

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: FITC-65244

Allgemeine Informationen

Katalog-Nr.:
FITC-65244

Größe:
100µg, 200 µg/ml

Wirt:
Maus

Isotyp:
IgG2a, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

MOPC-173

Anregungs-/Emissionsmaxima-
Wellenlängen:

495 nm / 524 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 0.2 ug FITC Plus Anti-Mouse IgG2a Isotype Control (FITC-65244, Clone: MOPC-173) (red), 0.2 ug FITC Plus Anti-Human CD8 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.