

Nur für Forschungszwecke

CD96 Monoklonaler Antikörper

Katalog-Nr.:FITC-65250



Allgemeine Informationen

Katalog-Nr.: FITC-65250	GenBank-Zugangsnummer: BC020749	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 µl/test	GeneID (NCBI): 10225	CloneNo.: NK92.39
Wirt: Maus	Vollständiger Name: CD96 molecule	Anregungs-/Emissionsmaxima-Wellenlängen: 495 nm / 524 nm
Isotyp: IgG1, kappa	Berechnete Masse: 402 aa, 46 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface co-stained with Anti-Human CD19 and 5 ul FITC Plus Anti-Human CD96 (FITC-65250, Clone: NK92.39), or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.