

Nur für Forschungszwecke

# CD33 Monoklonaler Antikörper

Katalog-Nr.: FITC-65272



## Allgemeine Informationen

<b>Katalog-Nr.:</b> FITC-65272	<b>GenBank-Zugangsnummer:</b> BC028152	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100tests , 5 µl/test	<b>GeneID (NCBI):</b> 945	<b>CloneNo.:</b> WM53
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> CD33 molecule	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 495 nm / 524 nm
<b>Isotyp:</b> IgG1, kappa	<b>Berechnete Masse:</b> 364 aa, 40 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

CD33, also referred to as Siglec-3, is a 67-kDa glycosylated transmembrane protein of sialic acid-binding immunoglobulin-like lectin (siglec) family. CD33 is expressed on monocytes, myeloid progenitors, granulocytes, mast cells, some T cells. It may mediate cell-to-cell adhesion and act as a receptor that inhibits the proliferation of normal and leukemic myeloid cells.

## Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Lagerungspuffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul FITC Plus Anti-Human CD33 (FITC-65272, Clone: WM53) (red) or unstained. Cells were not fixed. Monocytes were gated.