

Nur für Forschungszwecke

# DYKDDDDK tag Monoklonaler Antikörper

Katalog-Nr.: FITC-66008



## Allgemeine Informationen

Katalog-Nr.:  
FITC-66008

Größe:  
100ul, Konzentration: 1000 µg/ml von  
Nanodrop;

Wirt:  
Maus

Isotyp:  
IgG2b

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

8H6A10

Empfohlene Verdünnungen:

IF 1:50-1:500

Anregungs-/Emissionsmaxima-

Wellenlängen:

494 nm / 520 nm

## Anwendungen

Geprüfte Anwendungen:

IF

Getestete Reaktivität:

Rekombinantes Protein

Positivkontrollen:

IF : Transfizierte HEK-293-Zellen,

## Hintergrundinformationen

DYKDDDDK Tag ( Equivalent To FLAG Antibody From Sigma ) with the following sequence DYKDDDDK, is a hydrophilic tag for recombinant protein technology. Tags can be used as a tool to localize gene products in a variety of cell types, study proteins topology, and also help to identify and characterize new, low abundance or poorly immunogenic proteins. Due to its high hydrophilic character, the DYKDDDDK tag is likely to be located on the surface of a fusion protein, which enables the tag to be accessible for antibodies. DYKDDDDK Tag Antibody is generated against 1x DYKDDDDK tag (DYKDDDDK) and can recognize protein containing one or more DYKDDDDK tags, independently on N-terminal, C-terminal or internal regions of the target protein. Anti-FLAG is a registered trademark of Sigma-Aldrich Biotechnology.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Ethanol) fixed with Transfected HEK-293 cells using FITC-66008 (DDDDK tag antibody) at dilution of 1:50.