

Nur für Forschungszwecke

CD115 (c-fms) Monoklonaler Antikörper

Katalog-Nr.:PE-65095



Allgemeine Informationen

Katalog-Nr.:

PE-65095

Größe:

100ug, 0.2 mg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC036343

GeneID (NCBI):

12978

Vollständiger Name:

colony stimulating factor 1 receptor

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

AFS98

Anregungs-/Emissionsmaxima-Wellenlängen:

496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 0.13 ug PE Anti-Mouse CD115 (c-fms) (PE-65095, Clone: AFS98) (red) or 0.13 ug isotype control (black). Cells were not fixed.