

Nur für Forschungszwecke

CD24 Monoklonaler Antikörper

Katalog-Nr.:PE-65121



Allgemeine Informationen

Katalog-Nr.:
PE-65121

Größe:
100ug, 0.2 mg/ml

Wirt:
Ratte

Isotyp:
IgG2b

GenBank-Zugangsnummer:
BC075622

GeneID (NCBI):
12484

Vollständiger Name:
CD24a antigen

Reinigungsmethode:
Affinitätsreinigung

CloneNo.:
M1/69

Anregungs-/Emissionsmaxima-
Wellenlängen:
496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:
Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.06 ug PE Anti-Mouse CD24 (PE-65121, Clone: M1/69) (red) or 0.06 ug isotype control antibody (black). Cells were not fixed.