

Nur für Forschungszwecke

TER-119 Monoklonaler Antikörper

Katalog-Nr.: PE-65149



Allgemeine Informationen

Katalog-Nr.:

PE-65149

Größe:

100ug, 0.2 mg/ml

Wirt:

Ratte

Isotyp:

IgG2b

GenBank-Zugangsnummer:

GeneID (NCBI):

104231

Vollständiger Name:

lymphocyte antigen 76

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

TER-119

Anregungs-/Emissionsmaxima-
Wellenlängen:

496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

TER-119 is a mouse erythroid lineage-specific monoclonal antibody that reacts with erythroid cells at differentiation stages from early proerythroblast to mature erythrocyte, but not with cells showing typical erythroid blast-forming unit (BFU-E) and erythroid colony-forming unit (CFU-E) activities (PMID: 1975515; 10848813). TER-119 recognizes a 52-kDa molecule on erythrocyte membranes (PMID: 10848813). TER-119 antigen is a molecule associated with cell-surface glycophorin A but not with glycophorin A itself (PMID: 10848813).

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ C57BL/6 mouse bone marrow cells were surface stained with 0.50 ug PE Anti-Mouse TER-119 (PE-65149, Clone: TER-119) (red) or 0.50 ug isotype control antibody (blue). Cells were not fixed.