

Nur für Forschungszwecke

CD116 Monoklonaler Antikörper

Katalog-Nr.:PE-65160



Allgemeine Informationen

Katalog-Nr.:

PE-65160

Größe:

100tests, 10 µl/test

Wirt:

Maus

Isotyp:

Mouse IgG1 kappa

GenBank-Zugangsnummer:

BC002635

GeneID (NCBI):

1438

Vollständiger Name:

colony stimulating factor 2 receptor,
alpha, low-affinity (granulocyte-
macrophage)

Berechnete Masse:

46 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit R-
Phycoerythrin (PE) unter optimalen
Bedingungen konjugiert. Das
Konjugat wird mit Größenausschluss-
Chromatographie gereinigt.

CloneNo.:

4H1

**Anregungs-/Emissionsmaxima-
Wellenlängen:**

496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



0.1 mL human peripheral blood were surface stained with 10.00 ul/test PE Anti-Human CD116 (PE-65160, Clone:4H1) and then treated with red blood cell lysis buffer.