

Nur für Forschungszwecke

CD1a Monoklonaler Antikörper

Katalog-Nr.:PE-65172



Allgemeine Informationen

Katalog-Nr.:
PE-65172

Größe:
100tests , 20 µl/test

Wirt:
Maus

Isotyp:
Mouse IgG1

GenBank-Zugangsnummer:
BC031645

GeneID (NCBI):
909

Vollständiger Name:
CD1a molecule

Berechnete Masse:
327 aa, 37 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit R-Phycoerythrin (PE) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:
HI 149

Anregungs-/Emissionsmaxima-Wellenlängen:
496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern. Vor Licht schützen.
Lagerungspuffer:
PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ MOLT-4 cells were surface stained with 20.00 ul/test PE-Anti-Human CD1a antibody (PE-65172, Clone: HI149) (red) or PE-mouse IgG1 isotype control (blue). Cells were not fixed.