

Nur für Forschungszwecke

CD13 Monoklonaler Antikörper

Katalog-Nr.:PE-65186



Allgemeine Informationen

Katalog-Nr.:

PE-65186

Größe:

100tests , 5 µl/test

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC058928

GeneID (NCBI):

290

Vollständiger Name:

alanyl (membrane) aminopeptidase

Berechnete Masse:

110 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit R-Phycoerythrin (PE) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:

WM15

Anregungs-/Emissionsmaxima-Wellenlängen:

496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human, Nicht-Humane Primaten

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human whole blood granulocytes were surface stained with 5 ul PE Anti-Human CD13 (PE-65186, Clone:WM15) or Mouse IgG1 Isotype Control. Cells were not fixed.