

Nur für Forschungszwecke

IgG Isotype Control Monoklonaler Antikörper



Katalog-Nr.:PE-65210

Allgemeine Informationen

Katalog-Nr.:

PE-65210

Größe:

100ug, 0.2 mg/mL

Wirt:

Armenian Hamster

Isotyp:

IgG

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

PIP

Anregungs-/Emissionsmaxima-Wellenlängen:

496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Hintergrundinformationen

The PIP immunoglobulin is useful as an isotype-matched control. The PIP immunoglobulin is specific for Glutathione S-transferase and can be used as an isotype control for Armenian Hamster IgG antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.25 ug PE Armenian Hamster IgG Isotype Control (PE-65210, Clone: PIP) (blue) or 0.25 ug PE Anti-Mouse CD3ε (PE-65060, Clone: 145-2C11) (red). Cells were not fixed.