

Nur für Forschungszwecke

# IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: PE-65244

## Allgemeine Informationen

Katalog-Nr.:  
PE-65244

Größe:  
100ug, 0.2 mg/mL

Wirt:  
Maus

Isotyp:  
IgG2a, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

N/A

CloneNo.:

MOPC-173

Anregungs-/Emissionsmaxima-  
Wellenlängen:

496 nm, 565 nm / 578 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



$1 \times 10^6$  human PBMCs were surface stained with 0.2  $\mu$ g PE Anti-Mouse IgG2a Isotype Control (PE-65244, Clone:MOPC-173) (red), PE Anti-Human CD45RO (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.