

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-Marcksl1



Numéro de catalogue: 10034-1-AP

## Informations de base

Numéro de catalogue:  
10034-1-AP

Taille:  
150ul

Hôte:  
Lapin

Isotype:  
IgG

Numéro d'acquisition GenBank:  
BC006757

Identification du gène (NCBI):  
17357

Nom complet:  
MARCKS-like 1

MW calculé  
40 kDa

Méthode de purification:  
Purification par affinité contre  
l'antigène

## Applications

Applications testées:  
ELISA

Spécificité de l'espèce:  
Humain

## Informations générales

MacMARCKS is a member of MARCKS family of PKC substrate. It is widely used in the signal transduction studies as an indicator of PKC activation. The protein has a calculated molecular weight of 20 kDa and an apparent molecular weight of 42-45 kDa in different species on SDS PAGE. MacMARCKS functions in many aspects of cell physiology including integrin activation, cell spreading and phagocytosis. Its functions are thought to associate with its ability in binding and regulating actin and microtubule cytoskeleton. Current evidence indicated that this protein also involved controls membrane fluidity and molecular mobility of membrane proteins

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

