

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-TRAF4



Numéro de catalogue: 10083-2-AP

1 Publications

## Informations de base

Numéro de catalogue:

10083-2-AP

Taille:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG0121

Numéro d'acquisition GenBank:

BC001769

Identification du gène (NCBI):

9618

Nom complet:

TNF receptor-associated factor 4

MW calculé

54 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Demandes citées:

IHC

Spécificité de l'espèce:

Humain

Espèces citées:

Humain

## Informations générales

TRAF4 is the fourth member of the TNF receptor associated factor (TRAF) protein family. TRAF proteins are associated with and mediate signal transduction from members of the TNF receptor superfamily. TRAF1 and TRAF2 are two related proteins that form a heterodimeric complex that associates with the cytoplasmic domain of the TNF-receptor type 2. The third member TRAF3 associates with the cytoplasmic domain of CD40. TRAF4 has been shown to interact with neurotrophin receptor, p75 (NTR/NTSR1), and negatively regulate NTR induced cell death and NF-KAPPA B activation. This protein also binds to p47phox, a cytosolic regulatory factor included in a multi-protein complex known as NAD(P)H oxidase. TRAF4 thus, is thought to be involved in the oxidative activation of MAPK8/JNK.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Yu-Chan Chang	32188842	Cell Death Dis	IHC

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

