

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-KCNS3



Numéro de catalogue: 10092-1-AP

Informations de base

Numéro de catalogue:

10092-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG0137

Numéro d'acquisition GenBank:

BC004148

Identification du gène (NCBI):

3790

Nom complet:

potassium voltage-gated channel, delayed-rectifier, subfamily S, member 3

MW calculé

59 kDa

MW observés:

50 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules HeLa, tissu placentaire humain

Informations générales

KCNS3, also name as Kv9.3, belongs to S subfamily of potassium voltage-gated channel. It is highly expressed in pulmonary artery myocytes, placenta, and parvalbumin-containing GABA neurons in brain cortex. Its main functions are associated with the regulation of the resting membrane potential and the control of the shape and frequency of action potentials.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HeLa cells were subjected to SDS PAGE followed by western blot with 10092-1-AP (KCNS3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.