

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-Calpain 10



Numéro de catalogue: 10419-1-AP

Informations de base

Numéro de catalogue: 10419-1-AP	Numéro d'acquisition GenBank: BC004260	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 350 µg/ml by Nanodrop;	Identification du gène (NCBI): 11132	Dilutions recommandées: WB 1:500-1:1000
Hôte: Lapin	Nom complet: calpain 10	
Isotype: IgG	MW calculé 75 kDa	
Immunogen Catalog Number: AG0661	MW observés: 57 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules Jurkat,
Spécificité de l'espèce: Humain	

Informations générales

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Jurkat cells were subjected to SDS PAGE followed by western blot with 10419-1-AP (Calpain 10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.