

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PPAP2C



Numéro de catalogue: 10451-1-AP

Informations de base

Numéro de catalogue:

10451-1-AP

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG0660

Numéro d'acquisition GenBank:

BC002806

Identification du gène (NCBI):

8612

Nom complet:

phosphatidic acid phosphatase type 2C

MW calculé

32 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, rat, souris

Contrôles positifs:

WB : tissu hépatique de souris,

Informations générales

PPAP2C (phosphatidic acid phosphatase type 2C) is also named as PAP2C, LPP2, PAP-2c, PAP2-g. PPAP2C is a member of the phosphatidic acid phosphatase family of enzymes that regulate the dephosphorylation of lipid phosphates and the gene appears to be particularly highly expressed in breast carcinoma, dedifferentiated chondrosarcoma, fibromatosis, synovial sarcoma, and neurofibroma (PMID:19139135). It has 2 isoforms produced by alternative splicing and the full length protein has a glycosylation site.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse liver tissue were subjected to SDS PAGE followed by western blot with 10451-1-AP (PPAP2C Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.