

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-PTGES



Numéro de catalogue: 10596-1-AP

## Informations de base

Numéro de catalogue:  
10596-1-AP

Taille:  
150ul , Concentration: 133 µg/ml by  
Bradford method using BSA as the  
standard;

Hôte:  
Lapin

Isotype:  
IgG

Immunogen Catalog Number:  
AG0847

Numéro d'acquisition GenBank:  
BC008280

Identification du gène (NCBI):  
9536

Nom complet:  
prostaglandin E synthase

MW calculé  
17 kDa

Méthode de purification:  
Purification par affinité contre  
l'antigène

## Applications

Applications testées:  
ELISA

Spécificité de l'espèce:  
Humain, rat, souris

## Informations générales

Human prostaglandin E synthase (PTGES) is a member of a recently recognized protein superfamily consisting of membrane associated proteins involved in eicosanoid and glutathione metabolism (the MAPEG family) (PMID:10377395). It is also named as MGST1L1, MPGES1, PGES, PIG12 and belongs to the MAPEG family. This protein is expressed in membranes as an approximately 16-kD protein (PMID:10091672).

## Stockage

Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.  
Tampon de stockage:  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

