

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-GEMIN7



Numéro de catalogue: 10625-1-AP

## Informations de base

Numéro de catalogue:  
10625-1-AP

Taille:  
150ul , Concentration: 750 µg/ml by  
Nanodrop;

Hôte:  
Lapin

Isotype:  
IgG

Immunogen Catalog Number:  
AG0976

Numéro d'acquisition GenBank:  
BC007793

Identification du gène (NCBI):  
79760

Nom complet:  
gem (nuclear organelle) associated  
protein 7

MW calculé  
15 kDa

Méthode de purification:

Purification par affinité contre  
l'antigène

## Applications

Applications testées:  
ELISA

Spécificité de l'espèce:  
Humain, rat, souris

## Informations générales

### Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 10625-1-AP (GEMIN7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.