

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TOR3A



Numéro de catalogue: 10647-1-AP

Phare

Informations de base

Numéro de catalogue:

10647-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG1012

Numéro d'acquisition GenBank:

BC007571

Identification du gène (NCBI):

64222

Nom complet:

torsin family 3, member A

MW calculé

46 kDa

MW observés:

39-42 kDa, 46-52 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:300-1:1000

IF 1:50-1:500

Applications

Applications testées:

IF, WB, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

WB : cellules HEK-293, cellules HepG2, cellules PC-3, cellules Raji

IF : cellules HEK-293,

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




WB result of TOR3A antibody (10647-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-TOR3A transfected HEK-293 cells.



Various lysates were subjected to SDS PAGE followed by western blot with 10647-1-AP (TOR3A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



HEK-293 cells were subjected to SDS PAGE followed by western blot with 10647-1-AP (TOR3A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using 10647-1-AP (TOR3A antibody) at dilution of 1:50 and Alexa Fluor 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).