

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-DUSP13



Numéro de catalogue: 10909-1-AP

Informations de base

Numéro de catalogue:
10909-1-AP

Taille:
150ul, Concentration: 600 µg/ml by
Nanodrop;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG1348

Numéro d'acquisition GenBank:
BC009778

Identification du gène (NCBI):
51207

Nom complet:
dual specificity phosphatase 13

MW calculé
22 kDa, 36 kDa

MW observés:
22 kDa, 36 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:500-1:1000
IP 0.5-4.0 ug for IP and 1:500-1:1000
for WB
IHC 1:50-1:500

Applications

Applications testées:
IHC, IP, WB, ELISA

Spécificité de l'espèce:
Humain, rat, souris

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (*) À défaut, 'le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6.0.

Contrôles positifs:

WB : tissu testiculaire de souris, tissu de muscle squelettique humain, tissu testiculaire de rat

IP : tissu testiculaire de souris,

IHC : tissu de cancer du côlon humain,

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




Various lysates were subjected to SDS PAGE followed by western blot with 10909-1-AP (DUSP13 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



IP result of anti-DUSP13 (IP:10909-1-AP, 4ug; Detection:10909-1-AP 1:500) with mouse testis tissue lysate 4000 ug.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 10909-1-AP (DUSP13 antibody) at dilution of 1:100 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 10909-1-AP (DUSP13 antibody) at dilution of 1:100 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).