

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-GDE1



Numéro de catalogue: 11272-2-AP

1 Publications

Informations de base

Numéro de catalogue:

11272-2-AP

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG1821

Numéro d'acquisition GenBank:

BC014981

Identification du gène (NCBI):

51573

Nom complet:

glycerophosphodiester phosphodiesterase 1

MW calculé

38 kDa

MW observés:

55-66 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

Informations générales

GDE1(glycerophosphodiester phosphodiesterase 1) is also named as MIR16 and belongs to the glycerophosphoryl diester phosphodiesterase family. It is an 43 kDa integral membrane glycoprotein, which remains associated with membranes after alkaline (pH 11) extraction and is sensitive to endoglycosidase H and most likely involved in the remodeling of membrane PIs and regulation of membrane phospholipid composition, and hence in intracellular membrane trafficking, as well as in the modulation of G proteinmediated signaling(PMID:12576545). The full length protein has two glycosylation sites.

Publications notables

Autrice	Pubmed ID	Journal	Application
Sasha M Milbeck	32602905	Invest Ophthalmol Vis Sci	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

