

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ABI1



Numéro de catalogue: 11516-1-AP

Informations de base

Numéro de catalogue:
11516-1-AP

Taille:
150ul , Concentration: 533 µg/ml by
Bradford method using BSA as the
standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG2077

Numéro d'acquisition GenBank:
BC024254

Identification du gène (NCBI):
10006

Nom complet:
abl-interactor 1

MW calculé
507 aa, 55 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Applications

Applications testées:
IF, ELISA

Spécificité de l'espèce:
Humain, rat, souris

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of HepG2 cells, using ABL1 antibody 11516-1-AP at 1:100 dilution and PE-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).